

Exclusive pH electrode for the Pen type pH meter

# pH ELECTRODE

Model : PE-12

ISO-9001, CE, IEC1010

[www.yalab.com.tw](http://www.yalab.com.tw)

02-2389-0101



PE-12



PE-12

## FEATURES

- \* Exclusive pH electrode used for the Pen type pH meter PH-222, PH-223.
- \* Measuring Range :  
1 to 13 pH ( typical 0 to 14 pH ).
- \* Measuring Temperature :  
5 to 60 °C ( 41 to 140 °F ).
- \* Electrode Structure :  
Combination type.
- \* Zero Potential for  $7 \pm 1$  pH.
- \* Body Material : Epoxy.
- \* Build the Temp. sensor ( thermistor ), used for ATC ( automatic Temp. compensation ).



**Lutron**

[www.YaLAB.com.tw](http://www.YaLAB.com.tw)

*The Art of Measurement*