

meter to computer COM port interface cable

RS232 INTERFACE CABLE

Model : UPCB-01+UPCB-02+UPCB-04

ISO-9001, CE, IEC1010

www.yalab.com.tw 02-2389-0101



UPCB-01



UPCB-02



UPCB-04

UPCB-01, UPCB-02

* Available for all the LUTRON instruments
(build in the RS232 interface)

* The computer interface cable,
one side is the ear phone plug
that connect to the meter's RS232 output,
the other side is the D9 plug
that to connect the computer's COM terminal.

*UPCB-01

the direct interface cable.

*UPCB-02

the isolated (photo couple) interface cable.

UPCB-04

* RS232 cable

(Double female D9 RS232 cable)

available for

FG20KG-RS232 or FG5000A-RS232.



www.YaLAB.com.tw

The Art of Measurement

UPCB-01+UPCB-02+UPCB-04+0708