


7210-00 毛髮式 溫溼度記錄器 石英型

Thermohygrograph Quartz type

# TEMPERATURE HUMIDITY

**SK SATO**  
SATO KEIRYOKI MFG. CO., LTD.

## SIGMA II THERMOHYGROGRAPH



### SPECIFICATIONS

Cat. No	: No. 7210-00
Model	: NSII-Q
Clock	: Quartz type 2 "C" size batteries Battery life: about 1 year in continuous use
Drum rotation	: 1 day, 7 day and 32 day changed by speed switch
Sensors	: Aged bimetal strip for temperature Human hair bundle for relative humidity
Measuring range	: -15 to +40°C for temperature 0 to 100% for R. humidity
Accuracy	: ±1°C for temperature ±3% for R. humidity
Chart graduation	: 1°C for temperature 1% for R. humidity
Recording pens	: Cartridge pens (violet)
Dimensions	: 336(W) x 148(D) x 295(H) mm
Weights	: Approx. 3.1kg
Accessories	: 7 day chart for 1 year use (55 sheets/pad), chart holder, 2 "C" size batteries

**SK SATO KEIRYOKI MFG. CO., LTD.**

銷售商：邁多科技股份有限公司 電話：02-2389-0101 [www.YaLAB.com.tw](http://www.YaLAB.com.tw)