

optional ATC Temp. probe for pH meter

TEMPERATURE PROBE

Model : TP-06, TP-07

ISO-9001, CE, IEC1010

www.yalab.com.tw 02-2389-0101



TP-06

The optional Temp. probe (ATC probe) for pH meter PH-206.

Transistor sensor.

Range : 0 to 60 °C.

Accuracy : ± 0.8 °C / 1.5 °F.

TP-07

The optional Temp. probe (ATC probe) for pH meter

PH-207, PH-207HA, PH-208, PH-221, YK-2005WA, WA-2017SD, PH-221E, BPH-231.

Thermistor sensor.

Range : 0 to 60 °C.

Accuracy : ± 0.8 °C / 1.5 °F.



LT Lutron

www.YaLAB.com.tw

The Art of Measurement